

CDK4 Mouse mAb


ABZ4156

Product Overview


Name	CDK4 Mouse mAb
Accession(Primary)	P11802
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P11802)

Gene Name	CDK4
About Gene Symbol	CDK4 (Cyclin Dependent Kinase 4) is a protein kinase and signaling regulator involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of CDK4 in HeLa, Lncap and F9 lysates using CDK4 antibody.



Western blot analysis of CDK4 (3F9) in HeLa, C6 and F9 lysates using CDK4 antibody.