


Cyclin H (4E11) Mouse mAb

ABZ4157


Product Overview

Name	Cyclin H (4E11) Mouse mAb
Accession(Primary)	P51946
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (P51946)

Gene Name CCNH**About Gene Symbol**

Western blot analysis of Cyclin H in Jurkat, K562 and A431 lysates using Cyclin H antibody.



Immunoprecipitation analysis of Cyclin H (4E11)
in HeLa lysates using Cyclin H antibody.