

EEF2/Elongation factor 2 (4B3) Mouse mAb


ABZ4159

Product Overview


Name	EEF2/Elongation factor 2 (4B3) Mouse mAb
Accession(Primary)	P13639
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 95 kDa; Observed MW: 95 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; ChIP-1:20

Target information (P13639)

Gene Name	EEF2
About Gene Symbol	EEF2 (Eukaryotic Translation Elongation Factor 2) is a metabolic enzyme involved in translational regulation. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of eEF2 in HL-60, Jurkat, SHSY5Y, U20S and HeLa lysates using eEF2 antibody.



Immunocytochemistry analysis of Elongation factor 2 (4B3) in HeLa using eEF2 antibody.