

Syntenin 1 (3D9) Mouse mAb


ABZ4160

Product Overview

Name	Syntenin 1 (3D9) Mouse mAb
Accession(Primary)	O00560
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 32 kDa; Observed MW: 32 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (O00560)

Gene Name	SDCBP
About Gene Symbol	SDCBP (Syndecan Binding Protein) is a signaling scaffold protein involved in receptor-mediated signaling. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Syntenin in HeLa lysates
using Syntenin antibody.