

FYN (4B8) Mouse mAb

ABZ4161


Product Overview

Name	FYN (4B8) Mouse mAb
Accession(Primary)	P06241
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 61 kDa; Observed MW: 61 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200


Target information (P06241)

Gene Name FYN

About Gene Symbol



Western blot analysis of Fyn in HeLa, 3T3, C6, COS7, CEM, Ramos and Jurkat lysates using Fyn antibody.



Immunocytochemistry analysis of FYN (4B8) in HeLa using Fyn antibody.