

c-Jun Mouse mAb


ABZ4162

Product Overview

Name	c-Jun Mouse mAb
Accession(Primary)	P05412
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 43,48 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P05412)

Gene Name	JUN
About Gene Symbol	JUN (Jun Proto-Oncogene, AP-1 Transcription Factor Subunit) is a proto-oncogene involved in gene transcription and regulation. It is an important research target in cancer research and infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of c-Jun in CHO-K1
transfected by PEGFPC1 c-Jun and CHO-K1 lysates
using c-Jun antibody.