


Ezrin (2C12) Mouse mAb

ABZ4163

Product Overview

Name	Ezrin (2C12) Mouse mAb
Accession(Primary)	P15311
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 69 kDa; Observed MW: 81 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P15311)

Gene Name EZR**About Gene Symbol**

Western blot analysis of Ezrin in MCF-7, 3T3, A549 and rat Brain lysates using Ezrin antibody.

