


CKMT1 (1A6) Mouse mAb

ABZ4164

Product Overview

Name	CKMT1 (1A6) Mouse mAb
Accession(Primary)	P12532
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P12532)

Gene Name CKMT1A**About Gene Symbol**

Western blot analysis of CKMT1 in mouse Colon, 293 and A431 lysates using CKMT1 antibody.

