


NUP98 (3B8) Mouse mAb

ABZ4165


Product Overview

Name	NUP98 (3B8) Mouse mAb
Accession(Primary)	P52948
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 198 kDa; Observed MW: 98 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (P52948)

Gene Name NUP98**About Gene Symbol**

Western blot analysis of NUP98 in HeLa, Jurkat and MCF-7 lysates using NUP98 antibody.



Immunoprecipitation analysis of NUP98 (3B8) in HeLa lysates using NUP98 antibody.