


FYB (2H9) Mouse mAb

ABZ4166

Product Overview

Name	FYB (2H9) Mouse mAb
Accession(Primary)	O15117
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 85 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (O15117)

Gene Name FYB1**About Gene Symbol**

Western blot analysis of FYB in 10, 15 and 20ug Jurkat lysates using FYB antibody.

