


Lck (3E5) Mouse mAb

ABZ4168

Product Overview

Name	Lck (3E5) Mouse mAb
Accession(Primary)	P06239
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 58 kDa; Observed MW: 58 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P06239)

Gene Name LCK**About Gene Symbol**

Western blot analysis of LCK in Jurkat and Ramos lysates using LCK antibody.

