

Pyruvate Dehydrogenase E2 (4A4) Mouse mAb


ABZ4169

Product Overview


Name	Pyruvate Dehydrogenase E2 (4A4) Mouse mAb
Accession(Primary)	P10515
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 69 kDa; Observed MW: 69 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20

Target information (P10515)


Gene Name	DLAT
About Gene Symbol	DLAT (Dihydrolipoamide S-Acetyltransferase) is a metabolic dehydrogenase involved in adaptive immunity. It is an important research target in fibrotic disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Pyruvate Dehydrogenase E2 in Jurkat, A549, U251, F9, Lncap and HeLa lysates using Pyruvate Dehydrogenase E2 antibody.



Immunocytochemistry analysis of Pyruvate Dehydrogenase E2 (4A4) in HeLa using Pyruvate Dehydrogenase E2 antibody.



Immunoprecipitation analysis of Pyruvate Dehydrogenase E2 (4A4) in HeLa lysates using Pyruvate Dehydrogenase E2 antibody.

