

TXNRD2 (3F2) Mouse mAb


ABZ4170

Product Overview

Name	TXNRD2 (3F2) Mouse mAb
Accession(Primary)	Q9NNW7
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 57 kDa; Observed MW: 57 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q9NNW7)

Gene Name	TXNRD2
About Gene Symbol	TXNRD2 (Thioredoxin Reductase 2) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of TXNRD2 in A431, U251, U2OS, Lncap and HeLa lysates using TXNRD2 antibody.