

CKMT2 (3F4) Mouse mAb

ABZ4171


Product Overview

Name	CKMT2 (3F4) Mouse mAb
Accession(Primary)	P17540
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Rat
Calculated MW	Calculated MW: 48 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P17540)

Gene Name CKMT2

About Gene Symbol



Western blot analysis of CKMT2 in rat Brain lysates using CKMT2 antibody.

