

HDAC3 (3G3) Mouse mAb


ABZ4172

Product Overview


Name	HDAC3 (3G3) Mouse mAb
Accession(Primary)	O15379
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 49 kDa; Observed MW: 49 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (O15379)


Gene Name	HDAC3
About Gene Symbol	HDAC3 (Histone Deacetylase 3) is a histone deacetylase involved in chromatin remodeling and epigenetics. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of HDAC3 in THP-1, LNCAP, MCF-7, PC-12, 293, Jurkat, 3T3 and Hela lysates using HDAC3 antibody.



Western blot analysis of HDAC3 (3G3) in Jurkat and Hela lysates using HDAC3 antibody.



Western blot analysis of HDAC3 (3G3) in HDAC3 knock down 293T, wildtype (WT), K562, Jurkat and Hela lysates using HDAC3 antibody.