

# Cytochrome P450 17A1/CYP17A1 Mouse mAb

ABZ4173


## Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | Cytochrome P450 17A1/CYP17A1 Mouse mAb     |
| <b>Accession(Primary)</b>   | P05093                                     |
| <b>Host Species</b>         | Mouse                                      |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG2b                                      |
| <b>Species Reactivity</b>   | Human                                      |
| <b>Calculated MW</b>        | Calculated MW: 57 kDa; Observed MW: 60 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000                            |

## Target information (P05093)

**Gene Name** CYP17A1

**About Gene Symbol**



Western blot analysis of Cytochrome P450 17A1 in HeLa, HepG2 and A549 lysates using Cytochrome P450 17A1 antibody.

