


CBR1 Mouse mAb

ABZ4174


Product Overview

Name	CBR1 Mouse mAb
Accession(Primary)	P16152
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 30 kDa; Observed MW: 30 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (P16152)

Gene Name CBR1**About Gene Symbol**

Western blot analysis of CBR1 in HeLa, A431 and MDA-MB-468 lysates using CBR1 antibody.



Immunocytochemistry analysis of CBR1 (2C9) in Hela using CBR1 antibody.