


DR5 Mouse mAb

ABZ4177


Product Overview

Name	DR5 Mouse mAb
Accession(Primary)	O14763
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 48 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (O14763)

Gene Name TNFRSF10B**About Gene Symbol**

Western blot analysis of DR5 in mouse kidney, mouse brain, HCT116 and HeLa lysates using DR5 antibody.



Immunocytochemistry analysis of TNFRSF10B in HeLa cells using TNFRSF10B antibody.