


ALDH2 Mouse mAb

ABZ4179

Product Overview

Name	ALDH2 Mouse mAb
Accession(Primary)	P05091
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat, Cat, Hamster
Calculated MW	Calculated MW: 56 kDa; Observed MW: 56 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P05091)

Gene Name ALDH2**About Gene Symbol**

Western blot analysis of ALDH2 in HeLa, 3T3, C6, COS7, CHO-K1 and A549 lysates using ALDH2 antibody.

