

IMP3 (1F1) Mouse mAb


ABZ4180

Product Overview

Name	IMP3 (1F1) Mouse mAb
Accession(Primary)	O00425
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 64 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (O00425)

Gene Name	IGF2BP3
About Gene Symbol	IGF2BP3 (Insulin Like Growth Factor 2 mRNA Binding Protein 3) is a cytokine and immune signaling protein involved in adaptive immunity and translational regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of IMP3 in MCF-7, SKVO3, NTERA-2, Jurkat, HeLa and K562 lysates using IMP3 antibody.