


CD31 (2D3) Mouse mAb


ABZ4181

Product Overview

Name	CD31 (2D3) Mouse mAb
Accession(Primary)	P16284
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 83 kDa; Observed MW: 130 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (P16284)

Gene Name PECAM1**About Gene Symbol**



Immunoprecipitation analysis of CD31 (2D3) in THP-1 lysates using CD31 antibody.