

SMYD4 (7C8) Mouse mAb


ABZ4182

Product Overview


Name	SMYD4 (7C8) Mouse mAb
Accession(Primary)	Q8IYR2
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 89 kDa; Observed MW: 89 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q8IYR2)

Gene Name	SMYD4
About Gene Symbol	SMYD4 (SET And MYND Domain Containing 4) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of SMYD4 in HeLa, MCF-7 and A431 lysates using SMYD4 antibody.



Western blot analysis of SMYD4 (7C8) in HeLa, MCF-7 and MDA-MB-468 lysates using SMYD4 antibody.