


RPA32 (3E7) Mouse mAb

ABZ4183


Product Overview

Name	RPA32 (3E7) Mouse mAb
Accession(Primary)	P15927
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 29 kDa; Observed MW: 32 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20


Target information (P15927)

Gene Name RPA2**About Gene Symbol**

Western blot analysis of RPA32/RPA2 in MCF-7, HeLa and HUVEC lysates using RPA32/RPA2 antibody.



Immunocytochemistry analysis of RPA32 (3E7) in HeLa using RPA32/RPA2 antibody.



Immunoprecipitation analysis of RPA32 (3E7) in HeLa lysates using RPA32/RPA2 antibody.