

Hsp90 alpha/beta (4B6) Mouse mAb

ABZ4185


Product Overview

Name	Hsp90 alpha/beta (4B6) Mouse mAb
Accession(Primary)	P07900
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 84 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P07900)

Gene Name HSP90AA1/HSP90AB1

About Gene Symbol



Western blot analysis of Hsp90 in HeLa, 3T3, C6, K562 and A549 lysates using Hsp90 antibody.

