

Hsp90 beta (3B9) Mouse mAb


ABZ4186

Product Overview


Name	Hsp90 beta (3B9) Mouse mAb
Accession(Primary)	P08238
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat, Cat, Hamster
Calculated MW	Calculated MW: 83 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20

Target information (P08238)


Gene Name	HSP90AB1
About Gene Symbol	HSP90AB1 (Heat Shock Protein 90 Alpha Family Class B Member 1) is a molecular chaperone involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Hsp90 beta in HeLa, 3T3, C6, CHO-K1 and COS7 lysates using Hsp90 beta antibody.



Immunocytochemistry analysis of Hsp90 beta (3B9) in HeLa using Hsp90 beta antibody.



Immunoprecipitation analysis of Hsp90 beta (3B9)
in HeLa lysates using Hsp9 beta antibody.