

Ku70 (7B1) Mouse mAb


ABZ4187

Product Overview


Name	Ku70 (7B1) Mouse mAb
Accession(Primary)	P12956
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Cat
Calculated MW	Calculated MW: 70 kDa; Observed MW: 67 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20; ChIP-1:20

Target information (P12956)


Gene Name	XRCC6
About Gene Symbol	XRCC6 (X-Ray Repair Cross Complementing 6) is an ATP-dependent DNA helicase involved in DNA damage response and repair and adaptive immunity. It is an important research target in autoimmune disease and is widely studied to elucidate cellular mechanisms a



Western blot analysis of Ku70 in HeLa, A549, COS7 and K562 lysates using Ku70 antibody.



Immunocytochemistry analysis of Ku70 (7B1) in HeLa using Ku70 antibody.



Immunoprecipitation analysis of Ku70 (7B1) in HeLa lysates using Ku7 antibody.

