

Ku80 Mouse mAb


ABZ4188

Product Overview


Name	Ku80 Mouse mAb
Accession(Primary)	P13010
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Cat
Calculated MW	Calculated MW: 83 kDa; Observed MW: 86 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20; ChIP-1:20

Target information (P13010)


Gene Name	XRCC5
About Gene Symbol	



Western blot analysis of Ku80 in COS7, Jurkat, HeLa, MCF-7 and A549 lysates using Ku80 antibody.



Immunofluorescence analysis of Ku80 (8H1) in HeLa using Ku80 antibody.



Immunoprecipitation analysis of Ku80 (8H1) in
Hela lysates using Ku8 antibody.