

DDX3 Mouse mAb


ABZ4190

Product Overview


Name	DDX3 Mouse mAb
Accession(Primary)	O00571
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 73 kDa; Observed MW: 73 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20; ChIP-1:20

Target information (O00571)


Gene Name	DDX3X
About Gene Symbol	DDX3X (DEAD-Box Helicase 3 X-Linked) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of DDX3 in HeLa, 3T3, C6, COS7, K562 and Jurkat lysates using DDX3 antibody.



Immunocytochemistry analysis of DDX3 (6G8) in HeLa using DDX3 antibody.



Immunoprecipitation analysis of DDX3 (6G8) in
Hela lysates using DDX3 antibody.