


## XRCC1 Mouse mAb

ABZ4191

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | XRCC1 Mouse mAb                            |
| <b>Accession(Primary)</b>   | P18887                                     |
| <b>Host Species</b>         | Mouse                                      |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG1                                       |
| <b>Species Reactivity</b>   | Human                                      |
| <b>Calculated MW</b>        | Calculated MW: 69 kDa; Observed MW: 82 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000                            |

### Target information (P18887)

**Gene Name** XRCC1**About Gene Symbol**

Western blot analysis of XRCC1 in Raji, K562, HeLa, Jurkat and THP-1 lysates using XRCC1 antibody.

