


GADD45A (6C3) Mouse mAb

ABZ4194

Product Overview

Name	GADD45A (6C3) Mouse mAb
Accession(Primary)	P24522
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 18 kDa; Observed MW: 22 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P24522)

Gene Name GADD45A**About Gene Symbol**

Western blot analysis of GADD45 alpha in CHO-K1 lysates and CHO-K1 transfected by GADD45 alpha His fusion protein lysates using GADD45 alpha antibody.

