

BAG1 (2C12) Mouse mAb


ABZ4195

Product Overview

Name	BAG1 (2C12) Mouse mAb
Accession(Primary)	Q99933
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 39 kDa; Observed MW: 52,46,33 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q99933)

Gene Name	BAG1
About Gene Symbol	BAG1 (BAG Cochaperone 1) is a molecular chaperone involved in apoptosis. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of BAG1 (2C12) in CHO-K1 , transfected with a Human pEGFPC1BAG1 construct (A) or transfected with a Human pEGFPC1 construct (B), using Bag1 antibody.