


Rabbit IgG (5B2) Mouse mAb


ABZ4196

Product Overview


Name	Rabbit IgG (5B2) Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Rabbit
Calculated MW	Calculated MW: 55,25 kDa; Observed MW: 55,25 kDa
Recommended Dilution	WB-1:500-1:1000; ELISA-1:10000



Western blot analysis of MEK2 in THP-1, Jurkat, Raji, 3T3 and HeLa lysates using MEK2 antibody, then incubated with 201069(Rabbit IgG mAb). Anti-mouse IgG, HRPlinked antibody was used for detecting.



Western blot analysis of Rabbit IgG (5B2) in decreasing concentrations of total rabbit IgG reduced and denatured in 2X SDS loading buffer with DTT using IgG (HeavyChain Specific) antibody.



Immunoprecipitation analysis of Rabbit IgG (5B2) in CHO-K1 lysates transfected with HA-tagged NGAL using HATag antibody (lane 2). Lane 1 contains rabbit IgG IP control lysates. All lanes are probed with HATag antibody as the primary antibody. Panel A incubated with secondary antibodies anti-rabbit whole IgG, HRPlinked antibody for detecting. Panel B incubated with anti-rabbit IgG (HeavyChain Specific) mAb and the Anti-mouse whole IgG, HRPlinked antibody was used for detecting.