

eIF4E (8C8) Mouse mAb


ABZ4199

Product Overview


Name	eIF4E (8C8) Mouse mAb
Accession(Primary)	P06730
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 25 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P06730)

Gene Name	EIF4E
About Gene Symbol	EIF4E (Eukaryotic Translation Initiation Factor 4E) is a protein kinase and signaling regulator involved in translational regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of eIF4E in MDA-MB-468, ZR751, MCF-7, T47D and CHO-K1 lysates using eIF4E antibody.



Western blot analysis of eIF4E (8C8) in Jurkat, Hela, MCF-7, K562, C6, 3T3, PC-12, HepG2, CHO-K1, Raji, Raw264.7 and Ramos lysates using eIF4E antibody.