

## CREB (1E9) Mouse mAb


ABZ4201

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | CREB (1E9) Mouse mAb                       |
| <b>Accession(Primary)</b>   | P16220                                     |
| <b>Host Species</b>         | Mouse                                      |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG1                                       |
| <b>Species Reactivity</b>   | Human, Rat                                 |
| <b>Calculated MW</b>        | Calculated MW: 35 kDa; Observed MW: 43 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000                            |

### Target information (P16220)

|                          |                                                                                                                                                                                                                                                                |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | CREB1                                                                                                                                                                                                                                                          |
| <b>About Gene Symbol</b> | CREB1 (CAMP Responsive Element Binding Protein 1) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular |



Western blot analysis of CREB in Jurkat, A549, MCF-7, C6 and HeLa lysates using CREB antibody.