

4E BP3 (3B6) Mouse mAb

ABZ4204


Product Overview

Name	4E BP3 (3B6) Mouse mAb
Accession(Primary)	O60516
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 11 kDa; Observed MW: 11 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (O60516)

Gene Name EIF4EBP3

About Gene Symbol



Western blot analysis of EIF4EBP3 in CHO-K1 lysates and CHO-K1 transfected by EGFPEIF4EBP3 fragment lysates using EIF4EBP3 antibody.

