


S100A10 (6F4) Mouse mAb

ABZ4205


Product Overview

Name	S100A10 (6F4) Mouse mAb
Accession(Primary)	P08207
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 11 kDa; Observed MW: 11 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20


Target information (P08207)

Gene Name S100a10**About Gene Symbol**

Western blot analysis of S100A10 in C6, 3T3 and HeLa lysates using S100A10 antibody.



Immunocytochemistry analysis of S100A10 (6F4) in HeLa using S100A10 antibody.



Immunoprecipitation analysis of S100A10 (6F4) in HeLa lysates using S1A1 antibody.