

Toll-Like Receptor 3 (2F10) Mouse mAb


ABZ4206

Product Overview

Name	Toll-Like Receptor 3 (2F10) Mouse mAb
Accession(Primary)	O15455
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 104 kDa; Observed MW: Refer to figures
Recommended Dilution	WB-1:500-1:1000

Target information (O15455)

Gene Name	TLR3
About Gene Symbol	



Western blot analysis of Toll-Like Receptor 7 in CHO-K1 lysates and CHO-K1 transfected by Toll-Like Receptor 3 lysates using Toll-Like Receptor 3 antibody.