

SSB (4B10) Mouse mAb


ABZ4208

Product Overview

Name	SSB (4B10) Mouse mAb
Accession(Primary)	P05455
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 47 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P05455)

Gene Name	SSB
About Gene Symbol	SSB (Small RNA Binding Exonuclease Protection Factor La) is a protein kinase and signaling regulator involved in adaptive immunity. It is an important research target in autoimmune disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of La Antigen in MCF-7, A549, Jurkat and Hela lysates using La Antigen antibody.