

Hsp90 alpha(6H7) Mouse mAb

ABZ4212


Product Overview

Name	Hsp90 alpha(6H7) Mouse mAb
Accession(Primary)	P07900
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 85 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P07900)

Gene Name HSP90AA1

About Gene Symbol



Western blot analysis of Hsp90 alpha in R Brain, Jurkat, CHO-K1, COS7, C6, 3T3 and HeLa lysates using Hsp90 alpha antibody.

