


Cytokeratin 18 (6A8) Mouse mAb

ABZ4214


Product Overview

Name	Cytokeratin 18 (6A8) Mouse mAb
Accession(Primary)	P05783
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 48 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (P05783)

Gene Name KRT18**About Gene Symbol**

Western blot analysis of Cytokeratin 18 in A549 and HeLa lysates using Cytokeratin 18 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast carcinoma, hepatocarcinoma, stomach cancer and colon cancer tissue using CK18 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.