


Synaptophysin (10C8) Mouse mAb

ABZ4215


Product Overview

Name	Synaptophysin (10C8) Mouse mAb
Accession(Primary)	P08247
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100


Target information (P08247)

Gene Name SYP**About Gene Symbol**

Western blot analysis of Synaptophysin (10C8) in rat brain lysates using Synaptophysin (10C8) antibody.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Synaptophysin (10C8) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Synaptophysin (10C8) in rat Brain, rat Hippocampus, Hela, 3T3 lysates using Synaptophysin (10C8) antibody