

Cytokeratin 19 Mouse mAb

ABZ4216


Product Overview

Name	Cytokeratin 19 Mouse mAb
Accession(Primary)	P08727
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 44 kDa; Observed MW: 44 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100


Target information (P08727)

Gene Name KRT19

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded staining of Human colon cancer tissue using Cytokeratin 19 (10G10) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Cytokeratin 19 (10G10)
in PC-3, SKOV-3 and MCF-7 lysates using
Cytokeratin 19 (10G10) antibody