

PARP1 (10C2) Mouse mAb

ABZ4220


Product Overview

Name	PARP1 (10C2) Mouse mAb
Accession(Primary)	P09874
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 113 kDa; Observed MW: 116 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100


Target information (P09874)

Gene Name PARP1


About Gene Symbol




Western blot analysis of PARP1 (10C2) in HeLa, MCF-7 and Jurkat lysates using PARP (10C2) antibody.



Immunohistochemistry analysis of paraffin-embedded Human Tonsil Tissue using Cleaved PARP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of PARP1 (10C2) in K562, Ramos, Jurkat lysates using PARP (1C2) antibody.



Western blot analysis of PARP1 (10C2) in 3T3, HeLa lysates using PARP1 (10C2) antibody

