

Hsp70 1A (1D9) Mouse mAb

ABZ4222


Product Overview


Name	Hsp70 1A (1D9) Mouse mAb
Accession(Primary)	P0DMV8
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (P0DMV8)


Gene Name HSPA1A

About Gene Symbol







Immunofluorescence analysis of Hsp70 1A (1D9) in HeLa using Hsp70 1A (1D9) antibody(left) ,and DAPI (blue).



Western blot analysis of Hsp70 1A (1D9) in HeLa, mouse Brain lysates using Hsp70 1A (1D9) antibody.



Western blot analysis of Hsp70 1A (1D9) in A549, K562, HeLa lysates using Hsp70 1A (1D9) antibody



Immunohistochemistry analysis of paraffin-embedded Human tonsil tissue using HSP70 (1D9) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.