


## HP1 gamma (6D8) Mouse mAb

ABZ4226


### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | HP1 gamma (6D8) Mouse mAb                  |
| <b>Accession(Primary)</b>   | Q13185                                     |
| <b>Host Species</b>         | Mouse                                      |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG1                                       |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 21 kDa; Observed MW: 24 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-1:50-1:100            |


### Target information (Q13185)

**Gene Name** CBX3**About Gene Symbol**


Immunohistochemistry analysis of paraffin-embedded Human Colon Carcinoma Tissue using HP1 gamma (6D8) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using HP1 gamma (6D8) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of HP1 gamma (6D8) in MCF-7, Raw264.7, HeLa lysates using HP1 gamma (6D8) antibody



Western blot analysis of HP1 gamma (6D8) in CHO-K1, rat brain, C6, 3T3, HeLa lysates using HP1 gamma (6D8) antibody.

