

Hsp90 (5D6) Mouse mAb

ABZ4230


Product Overview

Name	Hsp90 (5D6) Mouse mAb
Accession(Primary)	P07900
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 85 kDa; Observed MW: 95 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100


Target information (P07900)

Gene Name HSP90AA1


About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human Colon Carcinoma using Hsp90 alpha (5D6) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Hsp90 alpha (5D6) in K562, Ramos, Jurkat, rat Brain lysates using Hsp90 alpha (5D6) antibody



Western blot analysis of Hsp90 alpha (5D6) in HeLa, A549, HL6, U2OS, C6 lysates using Hsp90 alpha (5D6) antibody