

ICAM1 (1C7) Mouse mAb

ABZ4236


Product Overview

Name	ICAM1 (1C7) Mouse mAb
Accession(Primary)	P05362
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human
Calculated MW	Calculated MW: 58 kDa; Observed MW: 96 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100


Target information (P05362)

Gene Name ICAM1

About Gene Symbol



Western blot analysis of CD54(ICAM1) in Raji and SW480 lysates using CD54(ICAM1) antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsils using CD54/ICAM1 (1C7) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.