


ICAM1 (4F5) Mouse mAb

ABZ4237


Product Overview

Name	ICAM1 (4F5) Mouse mAb
Accession(Primary)	P05362
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 58 kDa; Observed MW: 96 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (P05362)

Gene Name ICAM1**About Gene Symbol**

Western blot analysis of CD54(ICAM1) in SW480 and Raji lysates using CD54(ICAM1) antibody.



Immunocytochemistry analysis of ICAM1 (4F5) in SW480 using CD54(ICAM1) antibody.