


Cytokeratin 18 (2B5) Mouse mAb

ABZ4238


Product Overview

Name	Cytokeratin 18 (2B5) Mouse mAb
Accession(Primary)	P05783
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human, Cat
Calculated MW	Calculated MW: 48 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (P05783)

Gene Name KRT18**About Gene Symbol**

Western blot analysis of Keratin 18 in HeLa, COS7 and 293T lysates using Keratin 18 antibody.



Immunocytochemistry analysis of Cytokeratin 18 (2B5) in HeLa using Keratin 18 antibody.