

Hexokinase I Mouse mAb


ABZ4239

Product Overview

Name	Hexokinase I Mouse mAb
Accession(Primary)	P19367
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 102 kDa; Observed MW: 102 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P19367)

Gene Name	HK1
About Gene Symbol	HK1 (Brain Form Hexokinase) is a protein kinase and signaling regulator involved in cellular metabolism. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Hexokinase 1 in MCF-7, Jurkat, CHO-K1, C2C12, C6, 3T3 and HeLa lysates using Hexokinase 1 antibody.